



Becton Dickinson (CME) U.K.
Limited
Kincraig Business Park
Kincraig Road
Blackpool
FY2 0PJ
t: 01253 206700
bd.com/en-uk

4th March 2019

Attn: Purchasing and clinical departments

RE: T34 - End of Sales ref.

Affected infusion pump models: CME T34 Code 100-100PSMLTR

Dear valued customer,

At Becton Dickinson (CME) UK Limited we strive to maintain and upgrade the products and services we offer.

We would like to inform you that, due to regulation driven changes, the T34 ambulatory syringe pump ref. 100-100PSMLTR has undergone manufacturing improvements and updates. Compliance with the latest revision of the IEC 60601-1-8 Standard is a mandatory requirement since 14 April 2018 and remains a pillar of our infusion devices' CE marking.

BD would like to advise you that the T34 ambulatory syringe pump ref. 100-100PSMLTR has been discontinued and is no longer available for sale and has now been replaced with the T34 3.1 Edition ambulatory syringe pump ref. 999-103EN. We would like to take this opportunity to reassure you that your current 2nd edition T34 is still compliant for use. This document is to provide you with an overview of the latest updates. (Please see table 1 for details).

As part of this action BD will update your order internally to detail the new pump code. A customer service representative will contact you directly to confirm acceptance of the amendment to complete and expedite orders in progress as soon as possible, therefore it is not necessary that you submit a new purchase order. Please be assured pricing will not be affected by this update. We apologise for any additional administrative work this creates.


Table 1.

Becton Dickinson (CME) U.K. Limited
Registered in England and Wales: 4236707
Registered Office: Unit 1 Kincraig Business Park, Kincraig Road, Blackpool, Lancashire. FY2 0PJ

Advancing the world of
health



Becton Dickinson (CME) U.K. Limited
 Kinraig Business Park
 Kinraig Road
 Blackpool
 FY2 0PJ
 t: 01253 206700
bd.com/en-uk

		Overview of updates		
		Current	New	
 <p style="text-align: right;"><i>T34 3rd Edition</i></p>	Codes	100-100PSMLTR	999-103EN	
	Design updates Versus 2 nd Edition pump	PCBA to update 3.1 Edition alarms		
		PCBA to update single fault failures		
		Update DFU & TSM for 3.1 Edition compliance		
		Update device label IPX rating IP22		
		Lock out period capability once bolus is enabled		
		Update battery connectors		
		Update to front label universal symbols replace text		
		Insert patient guide for homecare		
		Update quick Reference guide for clinicians		
Carry Bag pouch label				

For customers with active service contracts, support will be available under the terms and conditions of your service contract.

Your local BD representative is available to assist in identifying an appropriate strategy to meet your ongoing needs. If you have any further questions or enquiries, please contact Customer Support on 01253 206700 or email customersupport@cmemedical.co.uk.

BD thanks you for your support.

Best Regards,

Heather Cameron

Product Manager, Ambulatory IV Pumps & Sets UKI